

Automatically detect each component of the PQRST complex, including onset, amplitude, and intervals.


Assess heart function with the Electrocardiography Module for LabChart using real-time recordings or saved ECG traces.

Use one of the presets designed for species-specific ECG cycles and beat ranges (Human, Guinea pig, Pig, Dog, Mouse, Rat, or Rabbit), or customize your own. Utilize the beat averaging feature to reduce noise and artifacts from any number of heartbeats. Exclude atypical or noisy waveforms with the ECG Beat Classifier.

Save ECG analysis independently of raw data. Visualize data with tabular and graphical PQRST complex start, end, amplitude, and interval reports. Export QT/RR, QT/time, RR/time, or waterfall ECG graphs.

Applications

- Pharmacokinetics
- Drug Trials
- Isolated Heart Studies
- Hemodynamic Measurements
- Cardiac Physiology



Left: Use Chart View to record ECG data and calculate QT interval in real-time.

Below: Use the Split Bar to view two data blocks at once, live or offline.



ECG Analysis Module

ECG Beat Classification

Categorize beats according to activity and isoelectric noise. Easily visualize QRS complexes and RR interval variance. Remove artifacts caused by movement, electrical interference, or baseline drift. Exclude atypical beats, such as extrasystole and supraventricular arrhythmias.

ECG Averaging View

Select which beats you want to average with ECG Averaging View and easily visualize the mean PQRST trace before and after an experimental intervention. Automatically label complex components or manually adjust labels as desired. Generate tabular data and graphs for individual traces, or mean trace start, end, amplitude, or PQRST interval.

ECG Table View

Calculate and record the parameters of each averaged ECG beat. Export this table to LabChart Data Pad or other software.

- RR, PR, QT, and QRS Intervals
- QTc
- T, P, and R Amplitudes

Analysis Plots



Generate scatter plots for presentation and posthoc analysis. Display selected data in several graph types:

- QT/RR Plot (QT Interval vs. RR Interval)
- QT/Time Plot (ideal for pharmacokinetics)
- RR/Time Plot (indicates HR variance)
- Waterfall Plot (3D waveform analysis)


Expedient workflow

Navigate to your data point of interest from any screen, with linked Chart, Beat Classification, Table, and Averaging Views. Quickly identify edited data marked in red.

Below: Visualize your data in automatically generated Analysis Plots.



Above: ECG Beat Classifier categorizes beats based on noise and form factor parameters.



Above: ECG Averaging View displays the mean PQRST trace from this Rabbit ECG, highlighting each wave component.

The figure shows the ECG Analysis Table View for a Rabbit ECG. It displays a table of ECG parameters for 18 beats. The table has columns for Time (s), RR Interval (s), Heart Rate (BPM), P Duration (s), QRS Interval (s), QT Interval (s), QTc (s), JT Interval (s), P Amplitude (mV), Q Amplitude (mV), and R Amplitude (mV). The status bar indicates Time 74.574 s, Beats 242 - 251 (using 10 of 10).

Time (s)	RR Interval (s)	Heart Rate (BPM)	P Duration (s)	QRS Interval (s)	QT Interval (s)	QTc (s)	JT Interval (s)	P Amplitude (mV)	Q Amplitude (mV)	R Amplitude (mV)		
1	0.718	0.305	193.2	0.05691	0.03600	0.1268	0.2276	0.09082	0.1788	-0.3338	1.739	
2	1.966	0.3135	191.4	0.05728	0.02777	0.04900	0.1269	0.2267	0.07790	0.1731	-0.3006	1.724
3	3.224	0.3138	191.2	0.05895	0.02797	0.04900	0.1259	0.2247	0.07765	0.1744	-0.2947	1.739
4	4.482	0.3137	191.2	0.05813	0.02568	0.03500	0.1260	0.2250	0.09102	0.1869	-0.3219	1.731
5	5.738	0.3080	194.8	0.05774	0.02557	0.04900	0.1256	0.2263	0.07658	0.1563	-0.3019	1.697
6	6.996	0.3090	194.2	0.05714	0.02463	0.03600	0.1265	0.2275	0.09045	0.1631	-0.2944	1.743
7	8.254	0.3085	194.5	0.06076	0.03069	0.03600	0.1271	0.2289	0.09111	0.1625	-0.2988	1.729
8	9.432	0.3123	192.2	0.05952	0.03002	0.04900	0.1267	0.2267	0.07785	0.1725	-0.3163	1.752
9	10.683	0.3127	191.8	0.05566	0.02240	0.04900	0.1267	0.2265	0.07868	0.1756	-0.3319	1.732
10	11.928	0.3083	194.6	0.05627	0.02418	0.03600	0.1268	0.2284	0.09079	0.1762	-0.3362	1.752
11	13.139	0.3083	194.6	0.05665	0.02292	0.04900	0.1258	0.2267	0.07785	0.1803	-0.3141	1.716
12	14.394	0.3033	197.0	0.05449	0.02139	0.04900	0.1259	0.2286	0.07809	0.2031	-0.2775	1.707
13	15.599	0.3003	199.8	0.06121	0.02851	0.03600	0.1268	0.2313	0.09075	0.2000	-0.3031	1.766
14	16.794	0.2985	207.3	0.05938	0.02688	0.04900	0.1242	0.2308	0.07517	0.1881	-0.3312	1.766
15	17.946	0.2985	201.0	0.05973	0.02638	0.03700	0.1255	0.2297	0.08847	0.1637	-0.3056	1.790
16	19.167	0.3127	191.8	0.05797	0.02512	0.03700	0.1257	0.2248	0.08870	0.1881	-0.2963	1.736
17	20.413	0.3105	193.2	0.06289	0.02930	0.03600	0.1250	0.2243	0.08951	0.1628	-0.2947	1.730
18	21.656	0.3095	193.8	0.05761	0.02408	0.03700	0.1261	0.2267	0.08910	0.1750	-0.3075	1.748

Above: Export data from Table View or add it to LabChart Data Pad.

Ordering Information

The ECG Analysis Module for LabChart can be purchased individually as an Add-On for LabChart 8 (MLS360/8, Windows), or as part of LabChart Pro (MLS260/8).

LabChart Pro includes LabChart software and all LabChart Modules, providing powerful data acquisition and analysis capabilities (available for Windows or Mac).

PowerLab and LabChart are trademarks of ADInstruments Pty Ltd. All other trademarks are the property of their respective owners. Products supplied by ADInstruments are intended for use in research and teaching applications and environments only.



Visit adstruments.com or contact your local ADInstruments representative for more information

Australia | Brazil | Europe | India | Japan | China | Middle East | New Zealand | North America | Pakistan | South America | South East Asia | United Kingdom

adstruments.com

